



1

- Place SPEEDX on a plane, smooth and hard surface
- Mettre SPEEDX sur une surface plane, lisse et rigide
- Apoya SPEEDX sobre una superficie llana, rigida y lisa
- Appoggia SPEEDX su una superficie piana, rigida e liscia
- SPEEDX auf einer ebenen, festen und Glatten oberfläche ablegen
- ПОМЕЩАТЬ SPEEDX НА ПЛОСКИЙ, ГЛАДКИЙ И ЖЁСТКИЙ ПЛОЩАДИ




2

- Position EASYPRESS
- Prendre EASYPRESS
- Coloca EASYPRESS
- Posiziona EASYPRESS
- EASYPRESS positionieren
- БРАТЬ EASYPRESS СО СЛЕДОМ ВНИЗ И НАПРАИТЬ В НУЖНОМ НАПРАВЛЕНИИ




3

- Centre EASYPRESS
- Positionner EASYPRESS
- Centra EASYPRESS
- Centra EASYPRESS
- EASYPRESS Zentrieren
- ФИКСИРОВАТЬ EASYPRESS




4

- Press EASYPRESS firmly
- Écraser avec force EASYPRESS
- Aprieta con fuerza EASYPRESS
- Premi con forza EASYPRESS
- EASYPRESS Kräftig druecken
- НАЖАТЬ С СИЛОЙ EASYPRESS




5

- Keep pressing EASYPRESS and massage immediately for 30 – 45 seconds
- En maintenant EASYPRESS écrasé, masser immédiatement durant 30-45 secondes
- Manten apretado EASYPRESS y fricciona inmediatamente durante 30-45 segundos
- Mantieni premuto EASYPRESS e massaggia immediatamente per 30-45 secondi
- EASYPRESS gedrückt halten und durchschnittlich 30-45 sekunden massieren
- ПОДДЕРЖИВАЯ НАЖАТЫМ EASYPRESS, МАССИРОВАТЬ БЫСТРО В ТЕЧНИЕ 30 – 45 СЕКУНД



6

- Push the liquid back
- Chasser le liquide vers la partie sachet
- Haz volver el liquido atras
- Fai ritornare il liquido indietro
- Flüssigkeit zurücklaufen lassen
- ПРОГОНЯТЬ ЖИДКОСТЬ К СТОРОНЕ ПАКЕТИК



7

- Open SPEEDX and take the film
- Ouvrir SPEEDX et extraire la pellicule
- Abre SPEEDX y extrae la pelicula
- Apri SPEEDX ed estrai la pellicola
- SPEEDX öffnen und den film entnehmen
- ОТКРЫВАТЬ SPEEDX И ИЗВЛЕКАТЬ ПЛЁНКУ



8

- Wash the developed film with care
- Bien laver la pellicule développée
- Lava bien la pelicula revelada
- Lava bene la pellicola sviluppata
- Den entwickelten film gut spülen
- МЫТЬ ПРОЯВЛЕННУЮ ПЛЁНКУ